

# BioFX™ Chemiluminescent Ultra Sensitive HRP Membrane Substrate

Product insert

### **Product number:**

CHMM-0060-2C (60 mL) Custom Sizes Available

## **Components:**

- 1) Chemiluminescent Ultra Sensitive Membrane Luminol-Reagent A (CHMC)
- 2) Chemiluminescent Ultra Sensitive Membrane Reaction Buffer-Reagent B (CHMD)

### Product use:

BioFX™ Chemiluminescent Ultra Sensitive HRP Membrane Substrate should be used for the sensitive detection of peroxidase in membrane blotting applications. Pat. www.surmodics.com/patents

# **Product stability and storage:**

Product stability	18 months from the date of manufacture.
Storage	Store at 2-8°C in the original container. The product should be protected from light.

### **Recommendations for use:**

- 1) Prepare the substrate by mixing 1 mL of Chemiluminescent Ultra Sensitive Membrane Luminol-Reagent A (CHMC) with 2 mL of Chemiluminescent Ultra Sensitive Membrane Reaction Buffer-Reagent B (CHMD).
  - The 1:2 ratio can be used to make desired volumes of substrate.
- 2) Mix the solutions well. The mixture of the two components is stable for several hours at room temperature when protected from light.

## Additional considerations:

The substrate is very sensitive and great care must be taken to optimize the individual assay components (antibodies, conjugates, solid phase, etc.) to minimize background reactivity associated with nonspecific immunochemical reactivities.

For technical assistance, email <a href="ivdtechsupport@surmodics.com">ivdtechsupport@surmodics.com</a>

## Related products:

In-Solution Protein Stabilizers & Diluents:	
StabilZyme™ HRP Conjugate Stabilizer (SZ02)	
StabilZyme™ SELECT Stabilizer (SZ03)	
StabilZyme™ NOBLE Stabilizer (SZ04)	
StabilZyme™ Protein Free Stabilizer (SZPF)	
Surmodics™ Assay Diluent (Protein-Free) (SM01)	
MatrixGuard™ Diluent (SM02)	



# BioFX™ Chemiluminescent Ultra Sensitive HRP Membrane Substrate

Product insert

Product Insert
Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
StabilBlot™ Protein-Free Blocker (SB01)
StabilBlot™ Milk Blocker – 10X (SB04)
StabilBlot™ Casein Blocker – 5X (SB05)
Substrates:
BioFX™ Chemiluminescent Ultra-Sensitive HRP Microwell Substrate (CHMI)
BioFX™ Chemiluminescent Sensitive Plus HRP Microwell and/or Membrane Substrate (LERI)
BioFX™ Chemiluminescent Super Sensitive HRP Microwell and/or Membrane Substrate (LUMI)
BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range One Component HRP Microwell Substrate (TMBX)
BioFX™ TMB Enhanced HRP Membrane Substrate (ESPM)
BioFX™ TMB One Component HRP Membrane Substrate (TMBM)

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.

19-May-2020 Product Code: CHMM Page 2
Doc # 7750, Revision 01