

Product number:

PNPS-0100-01 (100 mL)
PNPS-1000-01 (1000 mL)
Custom Sizes Available

Product use:

BioFX™ AP-Yellow One Component Microwell Substrate with Stabilizing Pellets contains *p*-nitrophenylphosphate (p-NPP) and is supplied as a one component solution for immunoassay applications. A soluble yellow product is obtained when this substrate is reacted with alkaline phosphatase. Pat. www.surmodics.com/patents

Product stability, storage and specifications:

Product stability	3 years from the date of manufacture.
Storage	2-8°C
Specification	This substrate is light sensitive and should be protected from direct ultraviolet sources.

Recommendations for use:

The product should be allowed to equilibrate to room temperature (25°C) prior to use. Unreacted substrate should be colorless to pale yellow in appearance and the removal of stabilizing pellets is not recommended.

- 1) Add 100-200 µL of substrate solution to each microwell.
- 2) A soluble yellow reaction product develops which can be read in the range of 405-420 nm.
 - For best results, sample absorbance values should be monitored and read before absorbance values exceed 2.0 OD units.
 - In endpoint assays, the substrate reaction can be stopped using an equal volume of AP Stop Reagent or 50 µL of 3 N NaOH per 200 µL of substrate.

Additional considerations:

Dilution of the substrate is not recommended. It is recommended that the antibodies or conjugates be diluted to reduce the intensity of the reaction.

The substrate is not recommended for membrane or immunohistochemical applications that require a precipitating reaction product.

For technical assistance, email ivdtechsupport@surmodics.com

Related products:

Aug-2020

Product Code: PNPS
Doc # 7765, Revision 01

Page 1

In-Solution Protein Stabilizers & Diluents:
StabilZyme™ AP Conjugate Stabilizer (SA01)
StabilZyme™ Protein Free AP Stabilizer (SAPF)
Surmodics™ Assay Diluent (Protein-Free) (SM01)
MatrixGuard™ Assay Diluent (SM02)
Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
Substrates:
BioFX™ AP Stop Solution (APSR)
BioFX™ Stabilizing Pellets for PNPS (STAB)

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.