

MPO, PR3 and BPI Antigens

Neutrophil granulocytes or polymorphonuclear leukocytes (PMN) are phagocytotic blood cells representing an important component of the mammalian innate immune system. Upon opsonization and subsequent phagocytosis by PMN, microorganisms are destroyed in the formed phagosome by a mechanism involving the production of reactive oxygen species and the secretion of proteolytic enzymes. Two important enzymes for this mechanism are myeloperoxidase (MPO) and proteinase 3 (PR3).

MPO is a heme containing enzyme composed of two identical subunits. Each subunit comprises a 64 kDa heavy and a 13 kDa light chain. MPO is stored in azurophilic granules, which fuse with and eventually release their content into the formed phagosome and the extracellular environment. MPO catalyzes the conversion of hydrogen peroxide to hypochlorous acid, a key step in the production of oxygen radicals leading to the oxidative degradation of microorganisms in the phagosomes. As an additional effect, hypochlorous acid produced in the extracellular environment is involved in the activation of proteases by inactivating protease inhibitors.

PR3 is a 29 kDa serine protease that is, like MPO, stored in azurophilic granules. In activated PMN, PR3 becomes secreted into both the phagosomes and the extracellular environment to destruct pathogens in the respective compartments.

Vasculitides are complex inflammatory autoimmune diseases of the blood vessels (vasculature) that appear to be associated with anti-neutrophil cytoplasmic antibodies (ANCA). Therefore, these diseases are also summarized as ANCA-associated vasculitides (AAV). In 1985, Woude *et al.* found an association between granulomatosis with polyangiitis (GPA; formerly known as Wegener's granulomatosis) and a subgroup of ANCA, which give rise to a characteristic cytoplasmic pattern (cANCA) in indirect immunofluorescence (IIF) using ethanol fixed granulocytes. Follow-up work by Niles *et al.* (1989) identified PR3 to be the antigen in cANCA. In 1988, Falk and Jennette described

a second subgroup of ANCA that appears to be associated with necrotizing and crescentic glomerulonephritis/renal vasculitis and give rise to a perinuclear or atypical cytoplasmic pattern (pANCA) in IIF using ethanol fixed granulocytes. MPO was identified as the antigen in pANCA.

However, later findings identified autoantibodies against other proteins that also yield in cANCA or pANCA in IIF making its interpretation more ambiguous. One example is the bactericidal/permeability-increasing protein (BPI), which is another component of azurophilic granules and is preferentially detected in a cANCA pattern (Zhao *et al.* 1995). Therefore, the international consensus statement published by Savige *et al.* (1999, 2003) recommended additional enzyme-linked immunosorbent assays (ELISA) to identify autoantibodies against MPO and PR3.

MPO, PR3, and BPI are purified from polymorphonuclear leukocytes of peripheral human blood.




Figure: Analyses of patient sera with confirmed cANCA and pANCA pattern in IIF by ELISA for the presence of autoantibodies against PR3 and MPO, respectively. As a control, sera from blood donors (BD) were included. The MPO positive cANCA serum is a diagnostic anomaly found in Churg-Strauss syndrome patients (R.-L. Humbel, personal communication).

References:

- Csernok and Moosig (2014) *Nat Rev Rheumatol.* 10 (8): 494-501
- Falk and Jennette (1988) *N Engl J Med.* 318 (25): 1651-1657
- Klebanoff (1968) *J Bacteriol.* 95 (6): 2131-2138
- Klebanoff *et al.* (2013) *J Leukoc Biol.* 93 (2): 185-198
- Korkmaz *et al.* (2010) *Pharmacol Rev.* 62 (4): 726-759
- Lionaki *et al.* (2012) *Arthritis Rheum.* 64 (10): 3452-3462
- Niles *et al.* (1989) *Blood.* 74 (6): 1888-1893
- Savige *et al.* (1999) *Am J Clin Pathol.* 111 (4): 507-513
- Savige *et al.* (2003) *Am J Clin Pathol.* 120 (3): 312-318
- Woude *et al.* (1985) *Lancet.* 1 (8426): 425-429
- Zhao *et al.* (1995) *Clin Exp Immunol.* 99 (1): 49-56

In some countries the use of certain antigens in diagnostic tests may be protected by patents. DIARECT is not responsible for the determination of these issues and suggests clarification prior to use.

Ordering Information

18500	Myeloperoxidase	0.1 mg
18501	(MPO; non recombinant)	1.0 mg
18600	Proteinase 3	0.1 mg
18601	(PR3; non recombinant)	1.0 mg
19200	BPI	0.1 mg
19201	(non recombinant)	1.0 mg

210416_Rev05

