

Product number:

ESPM-0100-01 (100 mL)
ESPM-1000-01 (1000 mL)
Custom Sizes Available

Product use:

BioFX™ TMB Enhanced One Component HRP Membrane Substrate is a ready-to-use solution containing 3,3',5,5'-tetramethylbenzidine in a mildly acidic buffer. The substrate reacts with peroxidase to form an insoluble permanent dark blue/purple color on the membrane or surface to which it is applied.

This enhanced formulation extends the dynamic range of detection, especially at higher concentrations of peroxidase where other substrates may exhibit loss of signal due to washout effects. Unreacted substrate should be faint yellow to pink in appearance.

Product stability, storage and specifications:

Product stability	Stable for 2 years from date of manufacture when stored at room temperature (25°C) or 4 years from date of manufacture when stored at 2-8°C.
Storage	2-25°C
Specification	This substrate is light sensitive and should be protected from direct UV sources. Over time, some flocculent material may appear; however, this does not affect product performance.

Recommendations for use:

The product should be allowed to equilibrate to room temperature (25°C) prior to use.

- 1) Cover the membrane surface with an ample amount of substrate and incubate for 5-20 minutes.
 - The substrate system will react with sites on the membrane containing peroxidase resulting in an insoluble permanent dark blue reaction product.
 - The reaction should be monitored and read before the background color becomes too intense which may result in diminished contrast between the positive and background staining.
- 2) Stop the reaction by rinsing the membrane with reagent quality water.

Additional considerations:

Dilution of the substrate is not recommended. It is recommended that the antibodies or conjugates be diluted to reduce the intensity of the reaction.

For technical assistance, email ivdtechsupport@surmodics.com



**BioFX™ TMB Enhanced One Component
HRP Membrane Substrate**

Product insert

Related products:

In-Solution Protein Stabilizers & Diluents:
StabilZyme™ HRP Conjugate Stabilizer (SZ02)
StabilZyme™ SELECT Stabilizer (SZ03)
StabilZyme™ NOBLE Stabilizer (SZ04)
StabilZyme™ Protein Free Stabilizer (SZPF)
Surmodics™ Assay Diluent (Protein-Free) (SM01)
MatrixGuard™ Diluent (SM02)
Blockers/Stabilizers:
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
StabilBlot™ Protein-Free Blocker (SB01)
StabilBlot™ Milk Blocker – 10X (SB04)
StabilBlot™ Casein Blocker – 5X (SB05)
Substrates:
BioFX™ Chemiluminescent Ultra-Sensitive HRP Microwell Substrate (CHMI)
BioFX™ Chemiluminescent Sensitive Plus HRP Microwell and/or Membrane Substrate (LERI)
BioFX™ Chemiluminescent Super Sensitive HRP Microwell and/or Membrane Substrate (LUMI)
BioFX™ TMB Super Sensitive One Component HRP Microwell Substrate (TMBS)
BioFX™ TMB One Component HRP Microwell Substrate (TMBW)
BioFX™ TMB Conductivity One Component HRP Microwell Substrate (TMBC)
BioFX™ TMB Slow Kinetic One Component HRP Microwell Substrate (TMSK)
BioFX™ TMB Extended Range One Component HRP Microwell Substrate (TMBX)
BioFX™ TMB One Component HRP Membrane Substrate (TMBM)

SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD, STABILBLOCK and STABILBLOT are trademarks of Surmodics, Inc. and/or its affiliates.

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.