

## ***Yersinia enterocolitica* (O:9) YopB**

### **Antigen Specification**

**Product Number: 43200**

**Description:**

*Yersinia enterocolitica* (O:9) YopB protein.  
Recombinant antigen for in vitro research and manufacturing use only.

**Immunological function:**

Binds IgG-, IgM, and IgA-type human antibodies.

**Origin:**

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

**Expression construct:**

Full-length cDNA coding for *Yersinia enterocolitica* (O:9) YopB fused to a deca-histidine purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

**Calculated molecular weight:**

43 kDa

**Calculated isoelectric point:**

pH 6.71

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

**Immunological tests/Functionality:**

Immunodot analyses with positive/negative samples.

**Recommended buffer/storage and handling conditions:**

Recommendations for storage buffer: neutral to slightly alkaline pH and 0.02% SDS.  
Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2016 - 2021 DIARECT GmbH