

Ribosomal Phosphoprotein P0

Antigen Specification

Product Number: 14100

Description:

Human ribosomal phosphoprotein P0 (belongs to a group of antigens also called acidic ribosomal phosphoproteins or ribosomal P-proteins).

Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for the human ribosomal P0 phosphoprotein fused to a hexa-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with i: anti-P0 autoantibody-positive sample; ii: monoclonal anti-His-tag antibody.

Calculated molecular weight:

35 kDa

Calculated isoelectric point:

pH 6.2

Immunological tests/Functionality:

Standard ELISA test (checkerboard analysis of positive/negative samples); immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH; due to purification workup under denaturing conditions presence of up to 0.02% SDS (or similar detergents) may be required for maintaining solubility. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Coating concentration:

0.3–0.7 µg/mL (depending on the type of ELISA plate and coating buffer). Suitable for labeling of functional groups.

Copyright 2002 - 2021 DIARECT GmbH