

PCNA humAb IgG

Antibody Specification

Product Number: 36600

Host: GammaPrim® mouse Clonality: monoclonal Isotype: IgG 1 / kappa

Description:

Human chimeric monoclonal antibody raised against recombinant Proliferating Cell Nuclear Antigen/Cyclin (PCNA). For in vitro research and manufacturing use only.

Immunogen/Epitope:

Human Proliferating Cell Nuclear Antigen/Cyclin (PCNA).

Immunological tests/Functionality:

Membrane-based immunoassay and ELISA with DIARECT's PCNA antigen.

Origin:

Produced by heterohybridoma cell line (constructed by fusing peripheral B lymphocytes from a GammaPrim® mouse immunized with PCNA to a hetero-hybridoma fusion cell line).

Concentration for assay use:

Concentrations for use as standard and/or positive control reagent in immunodiagnostic assays have to be individually determined by user.

Recommended buffer, storage and handling conditions:

Recommendations for storage buffer: neutral pH and 10% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2019 - 2020 DIARECT GmbH