

SLA/LP humAb IgG

Antibody Specification

Product Number: 36200

Host: GammaPrim[®] mouse

Clonality: monoclonal

Isotype: IgG 1 / kappa

Description:

Human chimeric monoclonal antibody raised against recombinant SLA/LP. For in vitro research and manufacturing use only.

Immunogen/Epitope:

Human SLA/LP.

Immunological function:

ELISA and membrane-based immunoassay with DIARECT's SLA/LP antigen.

Origin:

Produced by heterohybridoma cell line (constructed by fusing peripheral B lymphocytes from a GammaPrim[®] mouse immunized with SLA/LP to a hetero-hybridoma fusion cell line).

Concentration for assay use:

Concentrations for use as standard and/or positive control reagent in immunodiagnostic assays have to be individually determined by user.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral pH and 10% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2019 - 2021 DIARECT GmbH