

Bartonella henselae 26 kDa

Antigen Specification

Product Number: 45100

Description:

Bartonella henselae (strain Houston-1) 26 kDa protein. Recombinant antigen for research use or manufacturing only.

Origin:

Recombinant. Expressed in *E. coli* bacterial cells.

Expression construct:

Full-length cDNA coding for *Bartonella* henselae 26 kDa fused to a decahistidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

25 kDa

Calculated isoelectric point:

pH 7.27

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH, and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2017 - 2018 DIARECT AG