

***Borrelia afzelii* DbpA**

Antigen Specification

Product Number: 40900

Description:

Borrelia afzelii (strain pKo) decorin-binding protein A (DbpA).
Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG- and IgM-type human antibodies.

Origin:

Recombinant. Expressed in *E. coli* bacterial cells.

Expression construct:

Full-length cDNA coding for *Borrelia afzelii* DbpA fused to a hexa-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

17 kDa

Calculated isoelectric point:

pH 6.61

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with Lyme disease positive/negative samples; suitable for LTT (lymphocyte transformation test).

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2011 - 2021 DIARECT GmbH