



## BioFX™ Stabilizing Pellets for PNPS (pNPP)

Product insert

### Product number:

STAB-0020-01 (20 g)  
STAB-0050-01 (50 g)  
STAB-0100-01 (100 g)  
Custom Sizes Available

### Product use:

BioFX™ Stabilizing Pellets for PNPS (pNPP) are treated pellets for use in stabilizing phenyl phosphates.

### Product stability, storage and specifications:

Product stability	5 years from the date of manufacture.
Storage	Product should be stored at room temperature (25°C).

### Recommendations for use:

After dispensing BioFX product code PNPS, add the appropriate amount Stabilizing Pellets for PNPS (pNPP). Add 2 grams stabilizing pellets to 100 mL of PNPS substrate. Pellets should remain intact during measurement.

For technical assistance, email [ivdtechsupport@surmodics.com](mailto:ivdtechsupport@surmodics.com)

### Related products:

<b>In-Solution Protein Stabilizers &amp; Diluents:</b>
StabilZyme™ AP Conjugate Stabilizer (SA01)
StabilZyme™ Protein Free AP Stabilizer (SAPF)
Surmodics™ Assay Diluent (Protein-Free) (SM01)
MatrixGuard™ Assay Diluent (SM02)
<b>Blockers/Stabilizers:</b>
StabilGuard™ Immunoassay Stabilizer-BSA-Free (SG01)
StabilCoat™ Immunoassay Stabilizer (SC01)
StabilBlock™ Immunoassay Stabilizer (ST01)
<b>Substrates:</b>
BioFX™ AP Yellow One Component Microwell Substrate (pNPP)
BioFX™ AP Stop Reagent (APSR)

**SURMODICS, BIOFX, STABILZYME, MATRIXGUARD, STABILCOAT, STABILGUARD and STABILBLOCK are trademarks of Surmodics, Inc. and/or its affiliates.**

These products are for research and/or further manufacture and are not intended for use in humans or therapeutic purposes. Sales are without any seller's warranty or representation, expressed or implied, by usage or otherwise; no claims beyond replacement of unacceptable material or refund of purchase price shall be allowed. All claims must be made within 30 days following date of delivery.