

Ribosomal Phosphoprotein P1

Antigen Specification

Product Number: 14200

Description:

Human ribosomal phosphoprotein P1 (belongs to a group of antigens also called acidic ribosomal phosphoproteins or ribosomal P-proteins). Recombinant antigen for research use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for the human ribosomal P1 phosphoprotein fused to a hexa-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with i: anti-P1 autoantibody-positive sera; ii: monoclonal anti-hexa-His-tag antibody.

Calculated molecular weight:

12,336 Dalton

Calculated isoelectric point:

pH 4.75

Immunological tests/Functionality:

Standard ELISA test (checkerboard analysis of positive/negative sera panels); immunodot test with positive/negative sera panels.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH; due to purification workup under denaturing conditions presence of up to 0.02% SDS (or similar detergents) may be required for maintaining solubility. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Coating concentration:

0.3–0.7 µg/ml (depending on the type of ELISA plate and coating buffer). Suitable for labeling of functional groups.

Copyright 2002 – 2014 DIALECT AG