

## Borreliella afzelii DbpA

### **Antigen Specification**

**Product Number: 40900** 

#### **Description:**

Borreliella afzelii (strain pKo) decorinbinding protein A (DbpA). Recombinant antigen for research use or manufacturing only.

#### Immunological function:

Binds IgG- and IgM-type human antibodies.

#### Origin:

Recombinant. Expressed in *E. coli* bacterial cells.

#### **Expression construct:**

Full-length cDNA coding for *Borreliella afzelii* DbpA fused to a hexa-histidine purification tag.

#### **Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

## Calculated molecular weight:

17 kDa

## **Calculated isoelectric point:**

pH 6.61

#### **Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

#### Immunological tests/Functionality:

Immunodot analyses with Lyme disease positive/negative samples; suitable for LTT (lymphocyte transformation test).

# Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2011 - 2018 DIARECT AG