

***Borrelia burgdorferi* p58**

Antigen Specification

Product Number: 41600

Description:

Borrelia burgdorferi (strain B31) surface antigen p58.
Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human antibodies.

Origin:

Recombinant. Expressed in *E. coli* bacterial cells.

Expression construct:

Full-length cDNA coding for *Borrelia burgdorferi* p58 fused to a deca-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

60 kDa

Calculated isoelectric point:

pH 6.11

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with Lyme disease positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided.

Copyright 2014 - 2020 DIARECT GmbH