

***Borrelia burgdorferi* DbpB**

Antigen Specification

Product Number: 40600

Description:

Borrelia burgdorferi decorin-binding protein B (DbpB). Recombinant antigen for research use or manufacturing only.

Calculated molecular weight:

19,353 Dalton

Calculated isoelectric point:

pH 9.48

Immunological function:

Binds IgG- and IgM-type human antibodies.

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Origin:

Recombinant. Expressed in *E. coli* bacterial cells.

Immunological tests/Functionality:

Immunodot test with Lyme disease positive/negative plasma; suitable for LTT (lymphocyte transformation test).

Expression construct:

Full-length cDNA coding for *Borrelia burgdorferi* DbpB fused to a hexa-histidine purification tag.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with human plasma.

Copyright 2011 - 2014 DIARECT AG