

***Borrelia burgdorferi* p100**

Antigen Specification

Product Number: 40100

Description:

Borrelia burgdorferi (strain B31) p100 protein (p100/p83).
Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG- and IgM-type human antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for *Borrelia burgdorferi* p100 fused to a hexa-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

78 kDa

Calculated isoelectric point:

pH 5.25

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with Lyme disease positive/negative samples; suitable for LTT (lymphocyte transformation test) and ELISPOT (enzyme-linked immunospot).

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2011 - 2021 DIARECT GmbH