

Aln g 1.0101

Allergen Specification

Product Number: 50700

Description:

Alnus glutinosa Aln g 1.0101. Recombinant allergen for in vitro research and manufacturing use only.

Immunological function:

Binds IgE-type human antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for Aln g 1.0101 fused to a deca-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

19 kDa

Calculated isoelectric point:

pH 6.34

Molar extinction coefficient:

7450; A_{280} (1 mg/mL) = 0.398

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Storage conditions: -70°C or below.
Repeated freeze/thaw cycles should be avoided

Copyright 2015 - 2021 DIARECT GmbH