

# Borreliella burgdorferi p30

# **Antigen Specification**

**Product Number: 42600** 

## **Description:**

Borreliella burgdorferi (strain B31) p30 protein. Recombinant antigen for research use or manufacturing only.

## Immunological function:

Binds IgG- and IgM-type human antibodies.

## Origin:

Recombinant. Expressed in *E.coli* bacterial cells.

## **Expression construct:**

Full-length cDNA coding for *Borreliella* burgdorferi p30 fused to a hexahistidine purification tag.

#### **Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

## Calculated molecular weight:

29 kDa

## **Calculated isoelectric point:**

pH 5.34

#### **Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

### Immunological tests/Functionality:

Immunodot analyses with Lyme disease positive/negative samples.

# Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.

Storage conditions: -70°C or below.

Repeated freeze/thaw cycles should be avoided.

Copyright 2017 - 2018 DIARECT AG