

## ***Yersinia enterocolitica* (O:9) YopE**

### **Antigen Specification**

**Product Number: 43700**

**Description:**

*Yersinia enterocolitica* (O:9) YopE protein. Recombinant antigen for research use or manufacturing only.

**Immunological function:**

Binds IgG-, IgM, and IgA-type human antibodies.

**Origin:**

Recombinant. Expressed in *E.coli* bacterial cells.

**Expression construct:**

Full-length cDNA coding for *Yersinia enterocolitica* (O:9) YopE fused to a deca-histidine purification tag.

**Biochemical tests:**

SDS-PAGE (purity > 80%); Western blot with patient sample.

**Calculated molecular weight:**

24 kDa

**Calculated isoelectric point:**

pH 6.66

**Endotoxin content:**

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

**Recommended buffer/storage and handling conditions:**

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2016 - 2018 DIARECT AG