

Threonyl-tRNA Synthetase (PL-7)

Antigen Specification

Product Number: 15600

Description:

Human threonyl-tRNA synthetase (EC 6.1.1.3; PL-7).

Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-type human auto-antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for the human threonyl-tRNA synthetase (PL-7) fused to a hexa-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 90%); Western blot with i: anti-PL-7 autoantibody-positive sample; ii: monoclonal anti-His-tag antibody.

Calculated molecular weight:

83 kDa

Calculated isoelectric point:

pH 6.3

Immunological tests/Functionality:

Standard ELISA test (checkerboard analysis of positive/negative samples); immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent. Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Coating concentration:

0.3-0.8 µg/mL (depending on the type of ELISA plate and coating buffer). Suitable for labeling of functional groups.

Copyright 2003 - 2021 DIARECT GmbH