

***Yersinia enterocolitica* (O:9) YopB**

Antigen Specification

Product Number: 43200

Description:

Yersinia enterocolitica (O:9) YopB protein.
Recombinant antigen for in vitro research and manufacturing use only.

Immunological function:

Binds IgG-, IgM, and IgA-type human antibodies.

Origin:

Recombinant. Expressed by recombinant baculovirus (*Autographa californica* multiple nuclear polyhedrosis virus; AcMNPV) infection of *Spodoptera frugiperda* Sf9 insect cells.

Expression construct:

Full-length cDNA coding for *Yersinia enterocolitica* (O:9) YopB fused to a deca-histidine purification tag.

Biochemical tests:

SDS-PAGE (purity > 80%); Western blot with patient sample.

Calculated molecular weight:

43 kDa

Calculated isoelectric point:

pH 6.71

Endotoxin content:

Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.

Immunological tests/Functionality:

Immunodot analyses with positive/negative samples.

Recommended buffer/storage and handling conditions:

Recommendations for storage buffer: neutral to slightly alkaline pH and 0.02% SDS.

Storage conditions: -70°C or below. Repeated freeze/thaw cycles should be avoided.

Copyright 2016 - 2020 DIARECT GmbH